

ERMP1 Antibody

Catalog No: #40050



Package Size: #40050 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ERMP1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Affinity purified
Applications	ELISA WB
Species Reactivity	Human;Mouse;Rat
Immunogen Type	Recombinant Protein
Immunogen Description	Rabbit polyclonal ERMP1(1)antibody was raised against a recombinant humanERMP1 protein 204-365 aa(BC031630)
Conjugates	Unconjugated
Target Name	ERMP1
Concentration	2.8mg/ml
Formulation	PBS, pH 7.4 with 0.02% sodium azide.
Storage	This product is stable for several weeks at 4°C as an undiluted liquid. Dilute only prior to immediate use. For extended storage, aliquot contents and freeze at -20°C or below. Avoid cycles of freezing and thawing. Expiration date is one (1) year from date of receipt.

Application Details

ELISA: 1:20000-1:80000 WB: 1:200-1:1000

Images



Published Papers

el at., miR-9-5p attenuates ischemic stroke through targeting ERMP1-mediated endoplasmic reticulum stress. In Acta Histochem on 2019 Nov by Chi L, Jiao D,et al..PMID:31500865, , (2019)

[PMID:31500865](https://pubmed.ncbi.nlm.nih.gov/31500865/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.