

ENKD1 antibody

Catalog No: #40407



Package Size: #40407-1 50ul #40407-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	ENKD1 antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification
Applications	WB, IHC
Species Reactivity	Human;Mouse
Immunogen Description	Fusion protein of human ENKD1
Conjugates	Unconjugated
Target Name	ENKD1
Other Names	C16orf48; DAKV6410
Accession No.	NCBI Protein#:BC008284
Calculated MW	39 kDa
Concentration	1.4mg/ml
Formulation	pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol
Storage	Store at -20°C/1 year


Application Details

WB dilution:1:500-1:2000IHC dilution:1:50-1:200

Images



Gel: 8%SDS-PAGE, Lysate: 40 µg, Lane 1-3: Human cerebrum tissue, Human cerebella tissue, Mouse brain tissue lysates, Primary antibody:ENKD1 Antibody at dilution 1/1150, Secondary antibody: Goat anti rabbit IgG at 1/5000 dilution, Exposure time: 3 minutes



The image is immunohistochemistry of paraffin-embedded Human thyroid cancer tissue using (ENKbody) at dilution 1/80.

Background

The exact function of ENKD1 remains unknown.

Published Papers

Tianyu Zhang;Yixuan Wang;Xiaotong Nie;Xiangrong Chen;Yueyi Jin;Lulu Sun;Ruqian Yang;Jie Wang;Wenqing Xu;Ting Song;Wei Xie;Xiangfeng Chen;Chaojun Li;Jun Zhou;Sijin Wu;Yan Li;Tianliang Li et al., ENKD1 modulates innate immune responses through enhanced geranylgeranyl pyrophosphate synthase activity., , (2025)

[PMID:40048432](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.