XRCC4 Monoclonal Antibody

Catalog No: #40455

Package Size: #40455-1 50ul #40455-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


$\overline{}$		4.0	
	escri	nti	าท
$\boldsymbol{\nu}$	COUL	Pur	ווע

Product Name	XRCC4 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Clone No.	MR4516
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IP
Species Reactivity	Human
Specificity	The antibody detects endogenous XRCC4 proteins.
Conjugates	Unconjugated
Target Name	XRCC4
Other Names	XRCC4; DNA repair protein XRCC4; X-ray repair cross-complementing protein 4
Accession No.	Swiss-Prot: Q13426NCBI Gene ID: 7518
SDS-PAGE MW	38kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western blotting: 1:2000
Immunoprecipitation:1:200

Images



1) Input: Hela Cell Lysate 2) IP product: IP dilute 1:200

Western blot analysis of 1) Hela, 2) 293T, diluted at 1:3000.



Published Papers

el at., PHF6 promotes non-homologous end joining and G2 checkpoint recovery. In EMBO Rep on 2020 Jan 7 by Warmerdam DO, Alonso-de Vega I, et al..PMID:31782600, , (2019)

PMID:31782600

el at., Global non-covalent SUMO interaction networks reveal SUMO-dependent stabilization of the non-homologous end joining complex. In Cell Rep on 2021 Jan 26 by Romi ?

Gonzı ±η z-Prieto, Karolin Eifler-Olivi,et al..PMID:33503430, , (2021)

PMID:33503430

el at., Multi-Pathway DNA Double-Strand Break Repair Reporters Reveal Extensive Cross-Talk Between End-Joining, Single Strand Annealing, and Homologous Recombination.,, (2021)

PMID:

el at., Multi-pathway DNA-repair reporters reveal competition between end-joining, single-strand annealing and homologous recombination at Cas9-induced DNA double-strand breaks. In Nat Commun on 2022 Sep 8 by Bert van de Kooij, Alex Kruswick, et al..PMID:36075911, , (2022) PMID:36075911

Note: This product is for in vitro research use only and is not intended for use in humans or animals.