

Active Caspase-3 Monoclonal Antibody

Catalog No: #40491



Package Size: #40491-1 50ul #40491-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Active Caspase-3 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB IHC IF IP
Species Reactivity	Hu Ms Rt Chicken
Specificity	The antibody detects endogenous Active Caspase-3 protein.
Immunogen Type	Recombinant Protein
Immunogen Description	Recombinant Protein
Conjugates	Unconjugated
Target Name	Active Caspase-3
Accession No.	Swiss-Prot: P42574NCBI Gene ID: 836
SDS-PAGE MW	17kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:500-1000


Immunohistochemistry: 1:100-200

Immunofluorescence: IF 1:50-200

Images



Western blot analysis of 1) Hela, 2) 3T3, 3) Rat Brain Tissue using Active Caspase-3 Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using Active Caspase-3 Monoclonal Antibody.

Published Papers

JuBo Wang;DanFeng Dong;MaoDe Wang;Ke Gao et al., IDH1 Overexpression Induced Chemotherapy Resistance and IDH1 Mutation Enhanced Chemotherapy Sensitivity in Glioma Cells in Vitro and in Vivo, , (2014)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.