

GAPDH Monoclonal Antibody

Catalog No: #40493



Package Size: #40493-1 50ul #40493-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	GAPDH Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	The antibody was affinity-purified from antiserum by affinity-chromatography
Applications	WB
Species Reactivity	Hu Rt Ms Mk Dg C Hm Rb Pg Sh
Specificity	The antibody detects endogenous GAPDH protein.
Target Name	GAPDH
Other Names	GAPDH; GAPD; CDABP0047; OK/SW-cl.12; Glyceraldehyde-3-phosphate dehydrogenase; GAPDH; Peptidyl-cysteine S-nitrosylase GAPDH
Accession No.	Swiss-Prot: P04406NCBI Gene ID: 2597
SDS-PAGE MW	37kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:2000-5000

Images



Western blot analysis of Hela, diluted at 1:3000.

Published Papers

el at., Sorting nexin 27 interacts with Fzd7 and mediates Wnt signalling. In Biosci Rep on 2016 Jan 7 by Lixiang Sun, Xiuyong Hu et al.. PMID: 26744382, (2016)

PMID:26744382

et al., AMOTL2^Δ knockdown promotes the proliferation, migration and invasion of glioma by regulating EMT and Akt nuclear localization. In Oncol Rep on 2021 Jul by Xingjie Chen, Yalin Lu, et al.. PMID:34036399, (2021)

[PMID:34036399](#)

et al., A multi-omics approach based on 1H-NMR metabolomics combined with target protein analysis to reveal the mechanism of RIAISs on cervical carcinoma. In Aging (Albany NY) On 2022 Sep 27 by Chai Yanlan, Aibibai Jielili et al.. PMID:36170024, (2023)

[PMID:36170024](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.