

Bcl-w Polyclonal Antibody

Catalog No: #40641



Package Size: #40641-1 50ul #40641-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Bcl-w Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Bcl-w Polyclonal Antibody detects endogenous levels of Bcl-w protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Bcl-w.
Conjugates	Unconjugated
Target Name	Bcl-w
Other Names	BCL2L2; BCLW; KIAA0271; Bcl-2-like protein 2; Bcl2-L-2; Apoptosis regulator Bcl-W
Accession No.	Swiss-Prot: Q92843NCBI Gene ID: 599
SDS-PAGE MW	25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Bcl-w Polyclonal Antibody

Published Papers

et al., Treponema denticola Induces Neuronal Apoptosis by Promoting Amyloid- β Accumulation in Mice. In Pathogens on 2022 Oct 5 by Linrui Wu, Xinyi Su, et al..PMID: 36297207, , (2022)

[PMID:36297207](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.