

C/EBP β Polyclonal Antibody

Catalog No: #40657



Package Size: #40657-1 50ul #40657-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	C/EBP β Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	C/EBP β Polyclonal Antibody detects endogenous levels of C/EBP β protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human C/EBP β around the non-phosphorylation site of T235.
Conjugates	Unconjugated
Target Name	C/EBP β
Other Names	CEBPB; LAP; TCF5; PP9092; CCAAT/enhancer-binding protein beta; C/EBP beta; Liver activator protein; Nuclear factor NF-IL6; Transcription factor 5; TCF-5
Accession No.	Swiss-Prot: P17676NCBI Gene ID: 1051
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images

Western Blot analysis of 293T cells using C/EBP β Polyclonal Antibody

Published Papers

Yanlu Xiong;Jie Lei;Miaomiao Wen;Yongfu Ma;Jinbo Zhao;Yahui Tian;Zitong Wan;Xiaoyan Li;Jianfei Zhu;Wenchen Wang;Xiaohong Ji;Ying Sun;Jie Yang;Jiao Zhang;Shaowei Xin;Yang Liu;Lintao Jia;Yong Han;Tao Jiang et al., CENPF (+) cancer cells promote malignant progression of early-stage TP53 mutant lung adenocarcinoma., , (2025)

[PMID:40044674](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.