C1QBP Polyclonal Antibody

Catalog No: #40661

Package Size: #40661-1 50ul #40661-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

$\overline{}$		4.0	
	escri	nti	าท
$\boldsymbol{\nu}$	COUL	Pur	ווע

Concentration	1mg/ml
SDS-PAGE MW	32kd
Accession No.	Swiss-Prot: Q07021NCBI Gene ID: 708
	matrix protein p32; p33
	mitochondrial; GC1q-R protein; Glycoprotein gC1qBP; C1qBP; Hyaluronan-binding protein 1; Mitochondrial
Other Names	C1QBP; GC1QBP; HABP1; SF2P32; Complement component 1 Q subcomponent-binding protein;
Target Name	C1QBP
Conjugates	Unconjugated
Immunogen Description	Synthesized peptide derived from the Internal region of human C1QBP.
Specificity	C1QBP Polyclonal Antibody detects endogenous levels of C1QBP protein.
Species Reactivity	Human;Mouse;Rat
Applications	WB IHC ELISA
	immunogen.
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific
Clonality	Polyclonal
Host Species	Rabbit
Product Name	C1QBP Polyclonal Antibody

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using C1QBP Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.				