

# C1QBP Polyclonal Antibody

Catalog No: #40661

Package Size: #40661-1 50ul #40661-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                                                                                                                     |
|-----------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | C1QBP Polyclonal Antibody                                                                                                                                                                                           |
| Host Species          | Rabbit                                                                                                                                                                                                              |
| Clonality             | Polyclonal                                                                                                                                                                                                          |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                                                                               |
| Applications          | WB IHC ELISA                                                                                                                                                                                                        |
| Species Reactivity    | Human;Mouse;Rat                                                                                                                                                                                                     |
| Specificity           | C1QBP Polyclonal Antibody detects endogenous levels of C1QBP protein.                                                                                                                                               |
| Immunogen Description | Synthesized peptide derived from the Internal region of human C1QBP.                                                                                                                                                |
| Conjugates            | Unconjugated                                                                                                                                                                                                        |
| Target Name           | C1QBP                                                                                                                                                                                                               |
| Other Names           | C1QBP; GC1QBP; HABP1; SF2P32; Complement component 1 Q subcomponent-binding protein; mitochondrial; GC1q-R protein; Glycoprotein gC1qBP; C1qBP; Hyaluronan-binding protein 1; Mitochondrial matrix protein p32; p33 |
| Accession No.         | Swiss-Prot: Q07021NCBI Gene ID: 708                                                                                                                                                                                 |
| SDS-PAGE MW           | 32kd                                                                                                                                                                                                                |
| Concentration         | 1mg/ml                                                                                                                                                                                                              |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                                                                                                             |
| Storage               | Store at -20°C/1 year                                                                                                                                                                                               |

## Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of MCF7 cells using C1QBP Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.