

COL2A1 Polyclonal Antibody

Catalog No: #40772



Package Size: #40772-1 50ul #40772-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	COL2A1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	COL2A1 Polyclonal Antibody detects endogenous levels of COL2A1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human COL2A1.
Conjugates	Unconjugated
Target Name	COL2A1
Other Names	COL2A1; Collagen alpha-1(II) chain; Alpha-1 type II collagen
Accession No.	Swiss-Prot: P02458NCBI Gene ID: 1280
SDS-PAGE MW	140kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using COL2A1 Polyclonal Antibody



Western Blot analysis of K562 cells using COL2A1 Polyclonal Antibody

Published Papers

el at., 3D printed PGCL@PLA/10CSPL composite scaffolds loaded with fibronectin 1 for intervertebral disc degeneration treatment, In Biofabrication on 2024 Dec 13 by Weilin Zhang, Siyuan Chen,et al..PMID:39668784, , (2024)
[PMID:39668784](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.