

## Gα t1 Polyclonal Antibody

Catalog No: #40992



Package Size: #40992-1 50ul #40992-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Gα t1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Gα t1 Polyclonal Antibody detects endogenous levels of Gα t1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Gα t1.
Conjugates	Unconjugated
Target Name	Gα t1
Other Names	GNAT1; GNATR; Guanine nucleotide-binding protein G(t) subunit alpha-1; Transducin alpha-1 chain
Accession No.	Swiss-Prot: P11488NCBI Gene ID: 2779
SDS-PAGE MW	36kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


## Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/40000.

Not yet tested in other applications.

## Images



Western Blot analysis of COLO205 cells using Gα t1 Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.