

HDAC6 Polyclonal Antibody

Catalog No: #40999

Package Size: #40999-1 50ul #40999-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HDAC6 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse
Specificity	HDAC6 Polyclonal Antibody detects endogenous levels of HDAC6 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human HDAC6 around the non-phosphorylation site of S22.
Conjugates	Unconjugated
Target Name	HDAC6
Other Names	HDAC6; KIAA0901; JM21; Histone deacetylase 6; HD6
Accession No.	Swiss-Prot: Q9UBN7NCBI Gene ID: 10013
SDS-PAGE MW	131kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of HepG2 cells using HDAC6 Polyclonal Antibody

Published Papers

el at., Proteomic characterization of the colorectal cancer response to chemoradiation and targeted therapies reveals potential therapeutic strategiesIn
Cell Rep Med On2023 Dec 19byY Li, B Wang et al..PMID:38086380, , (2023)
[PMID:38086380](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.