

Ikaros 3 Polyclonal Antibody

Catalog No: #41055

Package Size: #41055-1 50ul #41055-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ikaros 3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Mouse
Specificity	Ikaros 3 Polyclonal Antibody detects endogenous levels of Ikaros 3 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ikaros 3.
Conjugates	Unconjugated
Target Name	Ikaros 3
Other Names	IKZF3; ZNFN1A3; Zinc finger protein Aiolos; Ikaros family zinc finger protein 3
Accession No.	Swiss-Prot: Q9UKT9NCBI Gene ID: 22806
SDS-PAGE MW	62kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COLO205 cells using Ikaros 3 Polyclonal Antibody

Published Papers

el at., Anti inflammatory and metabolic reprogramming effects of MENK produce antitumor response in CT26 tumor bearing mice.In J Leukoc Biol on 2020 Jan 29. by Tuo Y, Zhang Z, et al..PMID:31994797, , (2020)

[PMID:31994797](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.