

IL-5 Polyclonal Antibody

Catalog No: #41062

Package Size: #41062-1 50ul #41062-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-5 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA IF
Species Reactivity	Human;Mouse;Rat
Specificity	IL-5 Polyclonal Antibody detects endogenous levels of IL-5 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human IL-5.
Conjugates	Unconjugated
Target Name	IL-5
Other Names	IL5; Interleukin-5; IL-5; B-cell differentiation factor I; Eosinophil differentiation factor; T-cell replacing factor; TRF
Accession No.	Swiss-Prot: P05113NCBI Gene ID: 3567
SDS-PAGE MW	30kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

IF 1:50-200

ELISA: 1/5000.

Not yet tested in other applications.

Images



Western Blot analysis of HuvEc cells using IL-5 Polyclonal Antibody

Published Papers

el at., Mechanisms of regenerative preconditioning in pigs with subtotal hepatectomies. In Cir Cir on 2022 by Orue-Echebarria, Javier Vaquero,et al..PMID:35944117, , (2022)

[PMID:35944117](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.