

## Olfactory receptor 2M7 Polyclonal Antibody

Catalog No: #41269



Package Size: #41269-1 50ul #41269-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	Olfactory receptor 2M7 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IF ELISA
Species Reactivity	Hu
Specificity	Olfactory receptor 2M7 Polyclonal Antibody detects endogenous levels of Olfactory receptor 2M7 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the C-terminal region of human Olfactory receptor 2M7.
Target Name	Olfactory receptor 2M7
Other Names	OR2M7; Olfactory receptor 2M7; Olfactory receptor OR1-58
Accession No.	Swiss-Prot: Q8NG81NCBI Gene ID: 391196
SDS-PAGE MW	34kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details


Western Blot: 1/500 - 1/2000.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/10000.

Not yet tested in other applications.

## Images



Western Blot analysis of HepG2 cells using Olfactory receptor 2M7 Polyclonal Antibody

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.