

Paxillin Polyclonal Antibody

Catalog No: #41323



Package Size: #41323-1 50ul #41323-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Paxillin Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse;Rat;Monkey
Specificity	Paxillin Polyclonal Antibody detects endogenous levels of Paxillin protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from human Paxillin around the non-phosphorylation site of Y31.
Conjugates	Unconjugated
Target Name	Paxillin
Other Names	PXN; Paxillin
Accession No.	Swiss-Prot: P49023NCBI Gene ID: 5829
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using Paxillin Polyclonal Antibody

Published Papers

Xiaoyu Chen;Chengxia Shu;Wenqiang Li;Qiangqiang Hou;Guangmei Luo;Kexin Yang;Xiaoxing Wu et al., Discovery of a Novel Src Homology-2 Domain Containing Protein Tyrosine Phosphatase-2 (SHP2) and Cyclin-Dependent Kinase 4 (CDK4) Dual Inhibitor for the Treatment of Triple-Negative Breast Cancer, , (2022)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.