

PI 3-kinase p110 α Polyclonal Antibody

Catalog No: #41337

Package Size: #41337-1 50ul #41337-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description


Product Name	PI 3-kinase p110 α Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	IF;WB;IHC;ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	PI 3-kinase p110 α Polyclonal Antibody detects endogenous levels of PI 3-kinase p110 α protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human PI 3-kinase p110 α .
Conjugates	Unconjugated
Target Name	PI 3-kinase p110 α
Other Names	PIK3CA; Phosphatidylinositol 4; 5-bisphosphate 3-kinase catalytic subunit alpha isoform; PI3-kinase subunit alpha; PI3K-alpha; PI3Kalpha; PtdIns-3-kinase subunit alpha; Phosphatidylinositol 4,5-bisphosphate 3-kinase 110 kDa catalytic subunit
Accession No.	Swiss-Prot: P42336NCBI Gene ID: 5290
SDS-PAGE MW	110kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

IF 1:50-200 WB 1:500 - 1:2000. IHC 1:100 - 1:300. ELISA: 1:40000.

Images

Western Blot analysis of HELA cells using PI 3-kinase p110 α Polyclonal Antibody



Western Blot analysis of HELA cells using PI 3-kinase p110α Polyclonal Antibody

Published Papers

el at., EGFR-Phosphorylated Platelet Isoform of Phosphofructokinase 1 Promotes PI3K Activation. IN Mol Cell. On 2018 Apr 19 by Lee JH, Liu R et al.. PMID:29677490, , (2018)

[PMID:29677490](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.