

Ribosomal Protein L17 Polyclonal Antibody

Catalog No: #41401



Package Size: #41401-1 50ul #41401-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein L17 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Human;Mouse;Rat;Monkey
Specificity	Ribosomal Protein L17 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L17 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein L17.
Conjugates	Unconjugated
Target Name	Ribosomal Protein L17
Other Names	RPL17; 60S ribosomal protein L17; 60S ribosomal protein L23; PD-1
Accession No.	Swiss-Prot: P18621NCBI Gene ID: 6139
SDS-PAGE MW	24kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of COS7 cells using Ribosomal Protein L17 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.