

Ribosomal Protein L35 Polyclonal Antibody

Catalog No: #41407



Package Size: #41407-1 50ul #41407-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Ribosomal Protein L35 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Ribosomal Protein L35 Polyclonal Antibody detects endogenous levels of Ribosomal Protein L35 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Ribosomal Protein L35.
Conjugates	Unconjugated
Target Name	Ribosomal Protein L35
Other Names	RPL35; 60S ribosomal protein L35
Accession No.	Swiss-Prot: P42766NCBI Gene ID: 11224
SDS-PAGE MW	18kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

Immunohistochemistry: 1/100 - 1/300.

ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of Jurkat cells using Ribosomal Protein L35 Polyclonal Antibody

Note: This product is for in vitro research use only and is not intended for use in humans or animals.