

Siah-1/2 Polyclonal Antibody

Catalog No: #41434



Package Size: #41434-1 50ul #41434-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Siah-1/2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB IHC IF ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	Siah-1/2 Polyclonal Antibody detects endogenous levels of Siah-1/2 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the Internal region of human Siah-1/2.
Conjugates	Unconjugated
Target Name	Siah-1/2
Other Names	SIAH1; HUMSIAH; E3 ubiquitin-protein ligase SIAH1; Seven in absentia homolog 1; Siah-1; Siah-1a; SIAH2; E3 ubiquitin-protein ligase SIAH2; Seven in absentia homolog 2; Siah-2; hSiah2
Accession No.	Swiss-Prot: Q8IUQ4/O43255NCBI Gene ID: 6477/6478
SDS-PAGE MW	34kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.


ELISA: 1/10000.

Not yet tested in other applications.

Images



Western Blot analysis of HeLa cells using Siah-1/2 Polyclonal Antibody



Western Blot analysis of RAT-MUSCLE cells using Siah-1/2 Polyclonal Antibody

Published Papers

el at., MYPT1 Downregulation by LipopolysaccharideSIAH1/2 E3 LigaseUbiquitin-Proteasomal Degradation Contributes to Colonic Obstruction of Hirschsprung Disease.In Cell Mol Gastroenterol Hepatol on 2020 by Zhao W, Wang P et al..PMID:31759145, , (2020)
[PMID:31759145](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.