

## VASP Polyclonal Antibody

Catalog No: #41529



Package Size: #41529-1 50ul #41529-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

|                       |                                                                                                                       |
|-----------------------|-----------------------------------------------------------------------------------------------------------------------|
| Product Name          | VASP Polyclonal Antibody                                                                                              |
| Host Species          | Rabbit                                                                                                                |
| Clonality             | Polyclonal                                                                                                            |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Applications          | WB IHC IF ELISA                                                                                                       |
| Species Reactivity    | Human;Mouse;Rat;Monkey                                                                                                |
| Specificity           | VASP Polyclonal Antibody detects endogenous levels of VASP protein.                                                   |
| Immunogen Type        | peptide                                                                                                               |
| Immunogen Description | Synthesized peptide derived from human VASP around the non-phosphorylation site of S157.                              |
| Conjugates            | Unconjugated                                                                                                          |
| Target Name           | VASP                                                                                                                  |
| Other Names           | VASP; Vasodilator-stimulated phosphoprotein; VASP                                                                     |
| Accession No.         | Swiss-Prot: P50552NCBI Gene ID: 7408                                                                                  |
| SDS-PAGE MW           | 46+50kd                                                                                                               |
| Concentration         | 1mg/ml                                                                                                                |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| Storage               | Store at -20°C/1 year                                                                                                 |

## Application Details

Western Blot: 1/500 - 1/2000.


Immunohistochemistry: 1/100 - 1/300.

Immunofluorescence: 1/200 - 1/1000.

ELISA: 1/20000.

Not yet tested in other applications.

## Images



Western Blot analysis of COS7 cells using VASP Polyclonal Antibody

## Published Papers

el at., Identification of Novel Biomarkers With Diagnostic Value and Immune Infiltration in Burn Injury In Front Genet on 2022 Mar 22 by

Sitong Zhou, Kangchun Wang,et al..PMID:35391790, , (2022)

[PMID:35391790](#)

---

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.