

ALK-1 Polyclonal Antibody

Catalog No: #41591



Package Size: #41591-1 50ul #41591-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	ALK-1 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Mouse;Rat
Specificity	ALK-1 Polyclonal Antibody detects endogenous levels of ALK-1 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human ALK-1.
Conjugates	Unconjugated
Target Name	ALK-1
Other Names	ACVRL1; ACVRLK1; ALK1; Serine/threonine-protein kinase receptor R3; SKR3; Activin receptor-like kinase 1; ALK-1; TGF-B superfamily receptor type I; TSR-I
Accession No.	Swiss-Prot: P37023NCBI Gene ID: 94
SDS-PAGE MW	56kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using ALK-1 Polyclonal Antibody

Published Papers

el at., *Borrelia burgdorferi* Co-Localizing with Amyloid Markers in Alzheimer's Disease Brain Tissues. In *J Alzheimers Dis* on 2022 by Alireza G Senejani, Jasmin Maghsoudlou, et al.. PMID: 34897095, , (2022)

[PMID:34897095](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.