

IL-1RI Polyclonal Antibody

Catalog No: #41670

Package Size: #41670-1 50ul #41670-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	IL-1RI Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB,IHC-p,IF
Species Reactivity	Human;Mouse;Rat
Specificity	IL-1RI Polyclonal Antibody detects endogenous levels of IL-1RI protein.
Immunogen Description	Synthesized peptide derived from the N-terminal region of human IL-1RI.
Conjugates	Unconjugated
Target Name	IL-1RI
Other Names	IL1R1; IL1R; IL1RA; IL1RT1; Interleukin-1 receptor type 1; IL-1R-1; IL-1RT-1; IL-1RT1; CD121 antigen-like family member A; Interleukin-1 receptor alpha; IL-1R-alpha; Interleukin-1 receptor type I; p80; CD121a
Accession No.	Swiss-Prot: P14778NCBI Gene ID: 3554
SDS-PAGE MW	68kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of MCF7 cells using IL-1RI Polyclonal Antibody

Published Papers

el at., Functional heterogeneity of cancer-associated fibroblasts with distinct neoadjuvant immunotherapy plus chemotherapy response in esophageal squamous cell carcinoma. In Biomark Res on 2024 Sep 27 by Jun Jiang, Chao Xu, et al.. PMID:39334513, , (2024)

[PMID:39334513](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.