

Frizzled-10 Polyclonal Antibody

Catalog No: #41686

Package Size: #41686-1 50ul #41686-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Frizzled-10 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA
Species Reactivity	Human;Mouse
Specificity	Frizzled-10 Polyclonal Antibody detects endogenous levels of Frizzled-10 protein.
Immunogen Type	peptide
Immunogen Description	Synthesized peptide derived from the N-terminal region of human Frizzled-10.
Conjugates	Unconjugated
Target Name	Frizzled-10
Other Names	FZD10; Frizzled-10; Fz-10; hFz10; FzE7; CD350
Accession No.	Swiss-Prot: Q9ULW2NCBI Gene ID: 11211
SDS-PAGE MW	65kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

Images



Western Blot analysis of K562 cells using Frizzled-10 Polyclonal Antibody

Published Papers

et al., SEPT9_v2, frequently silenced by promoter hypermethylation, exerts anti-tumor functions through inactivation of Wnt/ β -catenin signaling pathway via miR92b-3p/FZD10 in nasopharyngeal carcinoma cells. In Clin Epigenetics on 2020 Mar 5; by Jiang Y, Liu L et al.. PMID:32138771, , (2020) [PMID:32138771](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.