

FOXP3 Polyclonal Antibody

Catalog No: #41822



Package Size: #41822-1 50ul #41822-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	FOXP3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB ELISA IHC IF
Species Reactivity	Human;Mouse;Rat
Specificity	FOXP3 Polyclonal Antibody detects endogenous levels of FOXP3 protein.
Immunogen Description	Synthesized peptide derived from the C-terminal region of human FOXP3.
Conjugates	Unconjugated
Target Name	FOXP3
Other Names	FOXP3; IPEX; JM2; Forkhead box protein P3; Scurfin
Accession No.	Swiss-Prot: Q9BZS1NCBI Gene ID: 50943
SDS-PAGE MW	47kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details


Western Blot: 1/500 - 1/2000.

ELISA: 1/20000.

Not yet tested in other applications.

IHC: 1/100-1/300

Images



Western Blot analysis of PC12 cells using FOXP3 Polyclonal Antibody.

Published Papers

et al., Effects of phenanthrene on oxidative stress and inflammation in lung and liver of female rats. In Environ Toxicol on 2020 Jan by Ma H, Wang H, et al..PMID: 31456356, , (2020)

[PMID:31456356](#)

et al., CD155 Cooperates with PD-1/PD-L1 to Promote Proliferation of Esophageal Squamous Cancer Cells via PI3K/Akt and MAPK Signaling Pathways. In Cancers (Basel) on 2022 Nov 15 by Xiyang Tan, Jie Yang, et al..PMID:36428703, , (2022)

[PMID:36428703](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.