

## TFF3 Monoclonal Antibody

Catalog No: #42028



Package Size: #42028-1 50ul #42028-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	TFF3 Monoclonal Antibody
Host Species	Mouse
Clonality	Monoclonal
Purification	protein G purified
Applications	IHC
Species Reactivity	Human
Specificity	specific for TFF3 native and denatured forms
Immunogen Type	protein
Immunogen Description	Recombinant Human Trefoil factor 3
Conjugates	Unconjugated
Target Name	TFF3
Other Names	Intestinal trefoil factor, hITF, Polypeptide P1.B, hP1.B, TFF3, ITF, TFI
Accession No.	Swiss-Prot#: Q07654
Concentration	1.0mg/mL
Formulation	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, PH 7.4
Storage	Store at -20°C


## Application Details

Immunohistochemistry: 1:20 - 1:200

## Images



Immunohistochemical analysis of paraffin-embedded human colon cancer using #42028 at dilution of 1:200.



Immunohistochemical analysis of paraffin-embedded human small intestine using #42028 at dilution of 1:200.

## Background

Involved in the maintenance and repair of the intestinal mucosa. Promotes the mobility of epithelial cells in healing processes (motogen).

## Published Papers

Chu Ling;Zhang Junbo;Meng Ruoyu;Zhuang Quan et al., Evaluating the role of IER3+ macrophages in the prognosis of liver fibrosis by bulk and single-cell transcriptional analyses, , (2024)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.