

MKI67/ki67 Antibody

Catalog No: #43095

Package Size: #43095-1 50ul #43095-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	MKI67/ki67 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	Antigen affinity purification.
Applications	IHC
Species Reactivity	Hu
Specificity	The antibody detects endogenous levels of total MKI67/ki67 protein.
Immunogen Type	peptide
Immunogen Description	Synthetic peptide of human MKI67/ki67
Target Name	MKI67/ki67
Other Names	KIA; MIB-; MIB-1; PPP1R105
Accession No.	Swiss-Prot#: P46013Gene ID: 4288
Concentration	0.8mg/ml
Formulation	Rabbit IgG in pH7.4 PBS, 0.05% NaN ₃ , 40% Glycerol.
Storage	Store at -20°C


Application Details

Immunohistochemistry: 1:25-1:100

Images



Immunohistochemical analysis of paraffin-embedded Human gastric cancer tissue using #43095 at dilution 1/20.



Immunohistochemical analysis of paraffin-embedded Human esophagus cancer tissue using #43095 at dilution 1/20.

Background

This gene encodes a nuclear protein that is associated with and may be necessary for cellular proliferation. Alternatively spliced transcript variants have been described. A related pseudogene exists on chromosome X.

Published Papers

et al., Long non-coding RNA PVT1 promotes tumor progression by regulating the miR-143/HK2 axis In gallbladder cancer. In Mol Cancer on 2019 Mar 2 by Chen J, Yu Y, et al.. PMID: 30825877, , (2019)

[PMID:30825877](#)

et al., Function of Cell-Cycle Regulators in Predicting Silent Pituitary Adenoma Progression Following Surgical Resection. In Oncol Lett on 2017 Dec by Sung Hyun Park, Ji Hwan Jang, et al.. PMID: 29344143, , (2017)

[PMID:29344143](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.