

HA-Tag Monoclonal Antibody(1B10)

Catalog No: #44005



Package Size: #44005-1 50ul #44005-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	HA-Tag Monoclonal Antibody(1B10)
Host Species	Mouse
Clonality	Monoclonal
Isotype	IgG1,Kappa
Applications	WB IP IF
Specificity	The antibody detects C-terminal, internal, and N-terminal HA-tag fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	HA-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Western blotting: 1:500-1:2000


Immunoprecipitation: 1:100

Immunofluorescence: 1:50-1:100


Images



IP antibody use: 5ug HA Mouse IgG1 per ml Lysate, WB 1: 5000 1) untransfected 293 cell lysate 2) transfected 293 cell lysate with HA-tag fusion protein 3) IP (untransfected 293+anti-HA mAb+Protein G agarose) 4) IP (transfected 293+ normal Mouse IgG+Prote



IF analysis of 293 cells transfected with a HA-tag protein, 1:2000 dilution (blue DAPI, red anti-HA)



Western blot analysis using HA tag mAb against recombinant protein. lane 1 :(100 ng); lane 2 :(50 ng); lane 3 :(25 ng); lane 4 :(10 ng); lane 5 :(2.5 ng); (Expected MW is 50 kDa)

Published Papers

el at., Identification of c-di-GMP/FleQ-Regulated New Target Genes, Including *cyaA*, Encoding Adenylate Cyclase, in *Pseudomonas putida*. In *mSystems* on 2021 May 11 by

Yujie Xiao, Haozhe Chen, et al.. PMID:33975969, (2021)

[PMID:33975969](https://pubmed.ncbi.nlm.nih.gov/33975969/)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.