

## V5-Tag Monoclonal Antibody(Q31)

Catalog No: #44024



Package Size: #44024-1 50ul #44024-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	V5-Tag Monoclonal Antibody(Q31)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q31
Isotype	IgG1
Applications	WB
Specificity	The antibody detects C-terminal, internal, and N-terminal V5-tag fusion proteins.
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Target Name	V5-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of V5 recombinant protein, (Q31) diluted at  
1) 1:5000  
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.