

Myc-Tag Monoclonal Antibody(5D6)

Catalog No: #44026



Package Size: #44026-1 50ul #44026-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com


Description

| | |
|-----------------------|---|
| Product Name | Myc-Tag Monoclonal Antibody(5D6) |
| Host Species | Mouse |
| Clonality | Monoclonal |
| Clone No. | 5D6 |
| Isotype | IgG1 |
| Applications | WB |
| Specificity | Antibody can detect C-Myc fusion proteins. |
| Immunogen Type | Peptide |
| Immunogen Description | Synthetic Peptide |
| Target Name | Myc-Tag |
| Concentration | 1mg/ml |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Store at -20°C/1 year |

Application Details

Western blotting: 1:500-1:2000

Images



Note: This product is for in vitro research use only and is not intended for use in humans or animals.