

# Myc-Tag Monoclonal Antibody(3E8)

Catalog No: #44029

Package Size: #44029-1 50ul #44029-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                         |
|-----------------------|-------------------------------------------------------------------------|
| Product Name          | Myc-Tag Monoclonal Antibody(3E8)                                        |
| Host Species          | Mouse                                                                   |
| Clonality             | Monoclonal                                                              |
| Clone No.             | 3.00E+08                                                                |
| Isotype               | IgG1                                                                    |
| Applications          | WB                                                                      |
| Specificity           | The antibody can detect C-Myc fusion proteins.                          |
| Immunogen Type        | Peptide                                                                 |
| Immunogen Description | Synthetic Peptide                                                       |
| Target Name           | Myc-Tag                                                                 |
| Concentration         | 1mg/ml                                                                  |
| Formulation           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage               | Store at -20°C/1 year                                                   |

## Application Details

Western blotting: 1:500-1:2000

## Images



Western blot analysis of c-Myc Recombinant protein, diluted at 1:2000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.