

Avi-Tag Monoclonal Antibody(6E3)

Catalog No: #44058



Package Size: #44058-1 50ul #44058-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	Avi-Tag Monoclonal Antibody(6E3)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	6.00E+03
Isotype	IgG1
Applications	WB
Specificity	The antibody detects Avi-Tag recombinant protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	Avi-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of Avi-Recombinant protein, diluted at
1) 1:5000
2) 1:10000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.