

RFP-Tag Monoclonal Antibody(Mix)

Catalog No: #44063

Package Size: #44063-1 50ul #44063-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	RFP-Tag Monoclonal Antibody(Mix)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Mix
Isotype	IgG1
Applications	WB
Specificity	The antibody detects RFP tag fusion proteins.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	RFP-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:500-1:2000

Images



Western blot analysis of recombinant RFP protein, diluted at 1:10,000.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.