

# ECFP-Tag Monoclonal Antibody(8F6)

Catalog No: #44067



Package Size: #44067-1 50ul #44067-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)

Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

Product Name	ECFP-Tag Monoclonal Antibody(8F6)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	8F6
Isotype	IgG1
Applications	WB
Specificity	The antibody detects endogenous ECFP-Tag protein.
Immunogen Type	Recombinant protein
Immunogen Description	Recombinant Protein
Target Name	ECFP-Tag
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

## Application Details

Western blotting: 1:2000-1:5000

## Images



Western blot analysis of Recombinant ECFP Protein using ECFP-Tag Monoclonal Antibody.

Note: This product is for in vitro research use only and is not intended for use in humans or animals.