TNF α Monoclonal antibody(Q34)

Catalog No: #44073

Package Size: #44073-1 50ul #44073-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

Product Name	TNF α Monoclonal antibody(Q34)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	Q34
Isotype	lgG1
Applications	WB IHC IF
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous TNF α protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	TNF-α
Other Names	TNF-α
Accession No.	P01375
Calculated MW	17,25kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

Western blotting: 1:1000-1:3000
Immunohistochemistry: 1:50-1:200


IF 1:50-200

Images




Western blot analysis of Recombinant Human TNF a Protein with TNF α Mouse mAb at 1:2,000.

The picture was kindly provided by our customer



Huazhong University of Science and Technology (HUST)



Immunohistochemical analysis of paraffin-embedded Human Tonsil Tissue using TNF aMouse mAb diluted at 1:50

Published Papers

el at., Screening of antii ζ Ahronic nonbacterial prostatitis activity of different extractions of the aerial part of Glycyrrhiza uralensis, and network pharmacology research. In Biomed Rep on 2021 Dec by Haifan Liu, Jie Cui,et al..PMID:34667596, , (2021)

PMID:34667596

el at., Screening of antii ζ Ahronic nonbacterial prostatitis activity of different extractions of the aerial part of Glycyrrhiza uralensis, and network pharmacology research. In Biomed Rep on 2021 Dec by Haifan Liu, Jie Cui,et al..PMID:34667596, , (2021)

PMID:34667596

Note: This product is for in vitro research use only and is not intended for use in humans or animals.