Bax Mouse Monoclonal Antibody(1C1)

Catalog No: #44150

Package Size: #44150-1 50ul #44150-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


Description

Product Name	Bax Mouse Monoclonal Antibody(1C1)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	1C1
Isotype	lgG1
Applications	IHC
Species Reactivity	Human;Mouse;Rat
Specificity	Bax Mouse Monoclonal Antibody(1C1) detects endogenous levels of BAX
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide
Conjugates	Unconjugated
Target Name	Bax
Other Names	BAX
Calculated MW	20kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:100-1:200


Images



Immunohistochemical analysis of paraffin-embedded Rat Testis Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Kidney Tissue using Bax Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Bax Mouse mAb diluted at 1:200

Published Papers

JuBo Wang;DanFeng Dong;MaoDe Wang;Ke Gao el at., IDH1 Overexpression Induced Chemotherapy Resistance and IDH1 Mutation Enhanced Chemotherapy Sensitivity in Glioma Cells in Vitro and in Vivo, , (2014)

PMID:

Note: This product is for in vitro research use only and is not intended for use in humans or animals.