Collagen I Mouse Monoclonal Antibody(6E11)

Catalog No: #44271

Package Size: #44271-1 50ul #44271-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com


Description

Product Name	Collagen I Mouse Monoclonal Antibody(6E11)
Host Species	Mouse
Clonality	Monoclonal
Clone No.	B'B 6E11
Isotype	lgG1
Applications	IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous Collagen I protein
Immunogen Type	Peptide
Immunogen Description	Synthetic Peptide of Collagen I
Conjugates	Unconjugated
Target Name	Collagen I
Other Names	Collagen alpha-1(I) chain (Alpha-1 type I collagen)
Calculated MW	139kd
Concentration	1mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year


Application Details

Immunohistochemistry: 1:50-1:200

Images



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Skin Carcinoma Tissue using Collagen I Mouse mAb diluted at 1:200.

Published Papers

el at., Combination therapy of chitosan, gynostemma, and motherwort alleviates the progression of experimental rat chronic renal failure by inhibiting STAT1 activation. In Oncotarget. On 2018 Mar 20 by Wenxia Bai, 1 Shudong Wang et al.. PMID: 29643988, , (2018)

PMID:29643988

Note: This product is for in vitro research use only and is not intended for use in humans or animals.