Product Datasheet

GPR41 Antibody

Catalog No: #44332

Package Size: #44332-1 50ul #44332-2 100ul



Orders: order@signalwayantibody.com Support: tech@signalwayantibody.com

Description

Product Name	GPR41 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific
	immunogen.
Applications	WB
Species Reactivity	Human
Specificity	GPR41 Antibody detects endogenous levels of total GPR41
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GPR41
Conjugates	Unconjugated
Target Name	GPR41
Other Names	FFAR3, G-protein coupled receptor 41, FFA3R, Gpcr41, Gpr42, Free fatty acid receptor 3, G protein-coupled
	receptor 41, GPR41
Accession No.	Swiss-Prot#: O14843NCBI Gene ID:2865
Calculated MW	38kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000

Published Papers

el at., Fecal microbiota transplantation and replenishment of short-chain fatty acids protect against chronic cerebral hypoperfusion-induced colonic dysfunction by . In J Neuroinflammation on 2022 Dec 26 by Shao-Hua Su, Yi-Fang Wu,et al..PMID:36567333, , (2022)

PMID:36567333

Note: This product is for in vitro research use only and is not intended for use in humans or animals.