

OR52R1 Antibody

Catalog No: #44647

Package Size: #44647-1 50ul #44647-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	OR52R1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	OR52R1 Antibody detects endogenous levels of total OR52R1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human OR52R1
Conjugates	Unconjugated
Target Name	OR52R1
Other Names	Olfactory receptor 52R1 Olfactory receptor OR11-22
Accession No.	Swiss-Prot#: Q8NGF1NCBI Gene ID:119695
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis OR52R1 using K562 whole cell lysates

Published Papers

el at., Validated Impacts of N6-Methyladenosine Methylated mRNAs on Apoptosis and Angiogenesis in Myocardial Infarction Based on MeRIP-Seq Analysis. In Front Mol Biosci on 2022 Jan 28 by Yingjie Zhang, Wenjie Hua,et al..PMID:35155564, , (2022)

[PMID:35155564](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.