

PLCE1 Antibody

Catalog No: #44842



Package Size: #44842-1 50ul #44842-2 100ul

Orders: order@signalwayantibody.com

Support: tech@signalwayantibody.com

Description

Product Name	PLCE1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	PLCE1 antibody detects endogenous levels of total PLCE1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human PLCE1
Conjugates	Unconjugated
Target Name	PLCE1
Other Names	1-phosphatidylinositol-4 antibody 5-bisphosphate phosphodiesterase epsilon-1 antibody Pancreas-enriched phospholipase C antibody Phosphoinositide phospholipase C-epsilon-1 antibody Phospholipase C-epsilon-1 antibody PLC-epsilon-1 antibody Plce1 antibody PLCE1_HUMAN antibody
Accession No.	Swiss-Prot#: Q9P212NCBI Gene ID:51196
Calculated MW	259 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: □ 1:500-1:2000

Published Papers

el at., Oncogene PLCE1 may be a diagnostic biomarker and prognostic biomarker by influencing cell cycle, proliferation, migration, and invasion ability in hepatocellular carcinoma cell lines. In J Cell Physiol on 2020 Feb 9. by Wang XK, Liao XW, et al.. PMID: 32037547, , (2020)

[PMID:32037547](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.