

KKLC1 Antibody

Catalog No: #45529

Package Size: #45529-1 50ul #45529-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	KKLC1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	KKLC1 Antibody detects endogenous levels of total KKLC1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human KKLC1
Conjugates	Unconjugated
Target Name	KKLC1
Other Names	CXorf61, KKLC1
Accession No.	Swiss-Prot#: Q5H943NCBI Gene ID:203413
Calculated MW	13 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis KKLC1 using HepG2 whole cell lysates

Published Papers

el at., Establishment of cancer/testis antigen profiling based on clinicopathological characteristics in resected pathological stage III non-small cell lung

cancer.In Cancer Manag Res.On 2018 Jul 16 by Jin S, Cao S et al..PMID: 30038519, , (2018)

[PMID:30038519](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.