

## IFNA7 Antibody

Catalog No: #45623

Package Size: #45623-1 50ul #45623-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)


## Description

|                       |                                                                                                                                                             |
|-----------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name          | IFNA7 Antibody                                                                                                                                              |
| Host Species          | Rabbit                                                                                                                                                      |
| Clonality             | Polyclonal                                                                                                                                                  |
| Purification          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                                       |
| Applications          | WB                                                                                                                                                          |
| Species Reactivity    | Human                                                                                                                                                       |
| Specificity           | IFNA7 Antibody detects endogenous levels of total IFNA7                                                                                                     |
| Immunogen Type        | Peptide                                                                                                                                                     |
| Immunogen Description | A synthesized peptide derived from human IFNA7                                                                                                              |
| Conjugates            | Unconjugated                                                                                                                                                |
| Target Name           | IFNA7                                                                                                                                                       |
| Other Names           | IFNA7, IFN-alpha 7, IFN-alpha-7, IFN-alpha-J1, IFN-alphaJ, IFNA-J, Interferon alpha-J, Interferon alpha-J1, Interferon alpha-7, LeIF J, Interferon, alpha 7 |
| Accession No.         | Swiss-Prot#: P01567NCBI Gene ID:3444                                                                                                                        |
| Calculated MW         | 22 kDa                                                                                                                                                      |
| Concentration         | 1.0mg/mL                                                                                                                                                    |
| Formulation           | Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.                                                          |
| Storage               | Store at -20°C                                                                                                                                              |

## Application Details

Western blotting: 1:1000-3000

## Images



Western blot analysis IFNA7 using HepG2 whole cell lysates

---

Note: This product is for in vitro research use only and is not intended for use in humans or animals.