

HECW1 Antibody

Catalog No: #45953



Package Size: #45953-1 50ul #45953-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	HECW1 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human
Specificity	HECW1 Antibody detects endogenous levels of total HECW1
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HECW1
Conjugates	Unconjugated
Target Name	HECW1
Other Names	KIAA0322, NEDL1
Accession No.	Swiss-Prot#: Q76N89NCBI Gene ID:23072
Calculated MW	180kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Published Papers

el at., The E3 ubiquitin ligase HECW1 targets thyroid transcription factor 1 (TTF1/NKX2. 1) for its degradation in the ubiquitin-proteasome system. In Cell Signal on 2019 Jun by Jia Liu, Su Dong,et al..PMID: 30849519, , (2019)

[PMID:30849519](#)

el at., Two distinct E3 ligases, SCFFBXL19 and HECW1, degrade thyroid transcription factor 1 in normal thyroid epithelial and follicular thyroid carcinoma cells, respectively. In FASEB J on 2019 Sep by Jia Liu, Su Dong, et al..PMID: 31238008, , (2019)

[PMID:31238008](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.