

## GNG5 Antibody

Catalog No: #46007



Package Size: #46007-1 50ul #46007-2 100ul

Orders: [order@signalwayantibody.com](mailto:order@signalwayantibody.com)Support: [tech@signalwayantibody.com](mailto:tech@signalwayantibody.com)

## Description

Product Name	GNG5 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	GNG5 Antibody detects endogenous levels of total GNG5
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human GNG5
Conjugates	Unconjugated
Target Name	GNG5
Other Names	Guanine nucleotide-binding protein G(I)/G(S)/G(O) subunit gamma-5,GNGT5
Accession No.	Swiss-Prot#: Q9Y3K8/P63218NCBI Gene ID:2787
Calculated MW	8 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

## Application Details

Western blotting: 1:1000-3000

Note: This product is for in vitro research use only and is not intended for use in humans or animals.