

HOXC13 Antibody

Catalog No: #46018



Package Size: #46018-1 50ul #46018-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	HOXC13 Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	HOXC13 Antibody detects endogenous levels of total HOXC13
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human HOXC13
Conjugates	Unconjugated
Target Name	HOXC13
Other Names	HOXC13, Homeo box C13, Homeobox protein Hox-3G, Homeobox protein Hox-C13, Homeo box 3G, HOX3G, ECTD9, Homeobox C13, HOX3, NUP98/HOXC13
Accession No.	Swiss-Prot#: P31276NCBI Gene ID:3229
Calculated MW	35 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis HOXC13 using HepG2 whole cell lysates

Published Papers

el at., The inconsistent regulation of HOXC13 on different keratins and the regulation mechanism on HOXC13 in cashmere goat (*Capra hircus*).In BMC Genomics. On 2018 Aug 23 by Wang S, Luo Z et al..PMID: 30139327, , (2018)
[PMID:30139327](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.