

LRP1B Antibody

Catalog No: #46033



Package Size: #46033-1 50ul #46033-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	LRP1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse
Specificity	LRP1B Antibody detects endogenous levels of total LRP1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human LRP1B
Conjugates	Unconjugated
Target Name	LRP1B
Other Names	LRP1B, LRPDIT, LRP-1B, LRP-deleted in tumors, LRP-DIT
Accession No.	Swiss-Prot#: Q9NZR2NCBI Gene ID:53353
Calculated MW	515kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis LRP1B using NIH-3T3 whole cell lysates

Published Papers

el at., Downregulation of Low-density lipoprotein receptor-related protein 1B (LRP1B) inhibits the progression of hepatocellular carcinoma cells by

activating the endoplasmic reticulum stress signaling pathway. In Bioengineered on 2022 Apr by Zili Zhen, Zhemin Shen, et al..PMID:35389768, , (2022)

[PMID:35389768](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.