

RASGEF1B Antibody

Catalog No: #46172



Package Size: #46172-1 50ul #46172-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com


Description

Product Name	RASGEF1B Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB
Species Reactivity	Human;Mouse;Rat
Specificity	RASGEF1B Antibody detects endogenous levels of total RASGEF1B
Immunogen Type	Peptide
Immunogen Description	A synthesized peptide derived from human RASGEF1B
Conjugates	Unconjugated
Target Name	RASGEF1B
Other Names	GPIG4BEB~Ras-GEF domain-containing family member 1B
Accession No.	Swiss-Prot#: Q0VAM2NCBI Gene ID:153020
Calculated MW	55 kDa
Concentration	1.0mg/mL
Formulation	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Storage	Store at -20°C

Application Details

Western blotting: 1:1000-3000

Images



Western blot analysis RASGEF1B using COLO205 whole cell lysates

Published Papers

Menglong Wang;Tongrui Wang;Yang Liu;Lurong Zhou;Yuanping Yin;Feng Gu et al., Identification and study of mood-related biomarkers and potential

molecular mechanisms in type 2 diabetes mellitus., (2025)

[PMID:39915429](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.