

Imp3 Polyclonal Antibody

Catalog No: #46748

Package Size: #46748-1 50ul #46748-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description


Product Name	Imp3 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Applications	WB IHC
Species Reactivity	Human;Mouse;Rat
Specificity	The antibody detects endogenous Imp3 protein
Immunogen Type	protein
Immunogen Description	Synthetic peptide from human protein at AA range: 135-184
Conjugates	Unconjugated
Target Name	Imp3
Other Names	IMP3 C15orf12 MRPS4
Accession No.	Swiss-Prot:Q9NV31NCBI Gene ID:55272
Calculated MW	25 kDa
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C

Application Details


Western blotting: 1B£B1500-1B£B110000

Immunohistochemistry: IHC: 1:100-1:300


Images



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200



Immunohistochemical analysis of paraffin-embedded human-lung-cancer, antibody was diluted at 1:200

Published Papers

el at., A novel tumour enhancer function of Insulin-like growth factor II mRNA-binding protein 3 in colorectal cancerInCell Death DisOn2023 Apr 6byDavide Di Fusco?1,?Antonio Di Grazia et al..PMID:37024466, , (2023)
[PMID:37024466](#)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.