

Transmembrane protein 200A Antibody

Catalog No: #48081

Package Size: #48081-1 50ul #48081-2 100ul

Orders: order@signalwayantibody.comSupport: tech@signalwayantibody.com

Description

Product Name	Transmembrane protein 200A Antibody
Host Species	Mouse
Clone No.	4-B3
Purification	Protein affinity purified
Applications	WB, ICC, IHC
Species Reactivity	Human;Mouse
Immunogen Description	Recombinant protein within Human Transmembrane protein 200A aa 59-359 / 491.
Conjugates	Unconjugated
Other Names	KIAA1913 antibody T200A_HUMAN antibody Tmem200a antibody Transmembrane protein 200A antibody TTMC antibody Two transmembrane C antibody
Accession No.	Swiss-Prot#:Q86VY9
Calculated MW	47 kDa
Formulation	1*PBS (pH7.4), 0.2% BSA, 40% Glycerol. Preservative: 0.05% Sodium Azide.
Storage	Store at -20°C

Application Details

WB: 1:500

IHC: 1:50

ICC: 1:100

Background

Transmembrane protein 200A(TMEM200A) is multi-pass membrane protein, which belongs to the TMEM200 family. It is expressed in cerebellum.

References

- Hui A B Y et al. Identification of a novel homozygous deletion region at 6q23.1 in medulloblastomas using high-resolution array comparative genomic hybridization analysis. Clin Cancer Res 11:4707-4716 (2005).
- Carninci P et al. The transcriptional landscape of the mammalian genome. Science 309:1559-1563 (2005).

Note: This product is for in vitro research use only and is not intended for use in humans or animals.