

PLC β2 Polyclonal Antibody

Catalog No: #48323



Package Size: #48323-1 50ul #48323-2 100ul

Orders: order@signalwayantibody.com
Support: tech@signalwayantibody.com


Description

Product Name	PLC β2 Polyclonal Antibody
Host Species	Rabbit
Clonality	Polyclonal
Isotype	IgG
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Applications	WB;IHC
Species Reactivity	Human;Mouse;Rat
Specificity	PLC β2 Polyclonal Antibody detects endogenous levels of PLC β2 protein.
Immunogen Description	Synthesized peptide derived from the Internal region of human PLC β2.
Conjugates	Unconjugated
Other Names	PLCB2; 1-phosphatidylinositol 4; 5-bisphosphate phosphodiesterase beta-2; Phosphoinositide phospholipase C-beta-2; Phospholipase C-beta-2; PLC-beta-2
Calculated MW	150kD
Concentration	1 mg/ml
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C/1 year

Application Details

WB 1:500-2000;IHC-p 1:50-300

Images



Immunohistochemical analysis of paraffin-embedded human Squamous cell carcinoma of lung.

Published Papers

el at., Salicin alleviates periodontitis via Tas2r143/gustducin signaling in fibroblasts. In Front Immunol on 2024 Mar 28 by Zhiying Zhang, Zhiyan Zhou,et al..PMID:38605968, , (2024)

Note: This product is for in vitro research use only and is not intended for use in humans or animals.